



**asma-1™**  
Electronic asthma monitor



**copd-6™**  
Electronic COPD screener



**lung monitor**  
Smoking cessation aid



**lung age**  
Lung function home monitors

	PEF	FEV <sub>1</sub>	FEV <sub>6</sub>	FEF 25-75%	Ratio	Lung Age	Memory	Personal Best	GOLD	PDF Reports
asma-1	●	●	○	○	○	○	600	●	○	○
asma-1 usb, s, bt	●	●	○	○	○	○	600	●	○	●
copd-6	○	●	●	○	●	●	○	○	●	○
copd-6 usb, s, bt	○	●	●	○	●	●	○	○	●	●
lung monitor	○	●	●	○	●	○	200	●	○	○
lung monitor usb, s, bt	○	●	●	●	●	○	200	●	○	●
lung age	○	●	○	○	○	●	○	○	○	○

● Available    ○ Not Available

usb - USB    bt - Bluetooth®    s - Serial

asma-1™ and copd-6™ are trademarks of Vitalograph Ltd.

The Bluetooth combination mark is a registered trademark and is wholly owned by the Bluetooth SIG.

# 4000 Series